

# NAV-ENABLED ABLATION CATHETERS



Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi




Select from these topics

- RHYTHMIA HDx™ Mapping System
- RHYTHMIA Catheter Portfolio
- INTELLAMAP ORION Mapping Catheter
- RHYTHMIA and INTELLANAV

## RHYTHMIA HDx™ Mapping System

**The RHYTHMIA HDx Mapping System's cutting-edge technology sets the standard of performance in high-definition mapping.**

RHYTHMIA's automated rapid acquisition of high-density, high-resolution maps provides unparalleled clarity so that you can efficiently identify the ablation target even in the most complex substrate.



Only the combination of advanced noise-filtering hardware, a novel software algorithm that can process unlimited data and a high-resolution mapping catheter working together can achieve true high-definition mapping.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- RHYTHMIA HDx™  
Mapping System
- RHYTHMIA  
Catheter Portfolio
- INTELLAMAP ORION  
Mapping Catheter
- RHYTHMIA and  
INTELLANAV

## RHYTHMIA™ Versatile Catheter Portfolio

The Boston Scientific nav-enabled catheter portfolio is a versatile suite of catheter tip technology combined with magnetic tracking for accuracy and efficiency. This offering is optimized for use with the high-resolution cardiac mapping system, RHYTHMIA HDx, giving you clarity for any complexity.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

RHYTHMIA HDx™  
Mapping System


RHYTHMIA  
Catheter Portfolio

INTELLAMAP ORION  
Mapping Catheter

RHYTHMIA and  
INTELLANAV

## INTELLAMAP ORION™ High-Resolution Mapping Catheter

The world's first variable deployment basket catheter (3-22mm)



**0.4 mm<sup>2</sup>**

*flat, printed* **ELECTRODES**

Maximize contact and minimize far field.  
See the clearest map possible.

**8x** *smaller*  
**ELECTRODES**<sup>1</sup>

Get sharper, better-quality signals for  
precise localization of arrhythmias.

**3x** *more*  
**ELECTRODES**

Collect more data more rapidly  
with 64 electrodes.

*Unique, Iridium-Oxide (Ir-Ox)  
coated electrode design for  
optimal noise performance.*





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi



Select from these topics

- RHYTHMIA HDx™ Mapping System
- RHYTHMIA Catheter Portfolio
- INTELLAMAP ORION Mapping Catheter
- RHYTHMIA and INTELLANAV

## RHYTHMIA™ and INTELLANAV™



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST



*Unparalleled Clarity  
Cool Performance  
Confident Navigation*



*Cool Performance  
Confident Navigation*



*Unparalleled Clarity  
Performance  
Accuracy*



*Familiarity  
Performance  
Accuracy*





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



## Select from these topics


- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - IO
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi 01
- IntellaNav 01
- IntellaNav MiFi XP
- IntellaNav ST

## Magnetic Sensor


Embedded magnetic sensors delivered tracking accuracy of better than 1mm<sup>1</sup>

Leverages hybrid tracking system: simultaneously creates impedance field maps while mapping

INTELLANAV  
MIFI™ XP EXAMPLE  
External View



Internal View



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - IO
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Enables Mapping Creation with RHYTHMIA™

Map Creation with INTELLANAV MIFI™ XP





Image courtesy of Matt Ostrom, MD, Torrance Memorial Hospital, Torrance, CA.

Map Creation with INTELLANAV™ OI





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi



Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - IO
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Proven BLAZER™ Mechanical Design



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - IO
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Built from the BLAZER™ Platform - *Now with Navigation*

### Predictable Handling

Designed for consistent handling – whether mapping or delivering therapy.

- Precise micro-movements
- Torqueability
- In-plane steering
- Tip stability

### Dependable Performance


Catheter can be relied upon to perform consistently throughout long cases.

- Curve retention

### Proven Results

One of the top-selling therapeutic catheters.

- Over 1 million catheters sold globally





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - OI
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Construction - Conventional Catheters

**INTELLANAV™ ST & INTELLANAV MIFI™ XP are most similar to the BLAZER™ II HTD / BLAZER II XP handling**

Catheter	Distal Shaft Material	Distal Shaft Durability	Pushability
BLAZER II	Single layer Pellethane™	Standard	Standard
BLAZER II HTD BLAZER II XP	Single layer Pebax™	Standard	Standard
BLAZER PRIME	Two layer Pebax	Fiber reinforced	Support mandrel
INTELLATIP MIFI XP	Two layer Pebax	Fiber reinforced	Standard
INTELLANAV ST INTELLANAV MIFI XP	Two layer Pebax	Fiber reinforced	Standard



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - OI
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Construction - Open-Irrigated

**INTELLANAV™ OI is most similar to BLAZER™ II Standard handling**  
**INTELLANAV MIFI™ OI is most similar to BLAZER II HTD handling**

Catheter	Distal Shaft Material	Distal Shaft Durability	Pushability
BLAZER II STD	Single layer Pellethane™	Standard	Standard
BLAZER II HTD	Single layer Pebax™	Standard	Standard
BLAZER PRIME	Two layer Pebax	Fiber reinforced	Support mandrel
BLAZER OI INTELLANAV OI	Single layer Pellethane	Standard	Standard
INTELLATIP MIFI OI INTELLANAV MIFI OI	Two layer Pebax	Fiber reinforced	Standard



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - OI
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## INTELLANAV™ Specification Summary

### INTELLANAV Ablation Catheters

	INTELLANAV ST	INTELLANAV MIFI XP	INTELLANAV OI	INTELLANAV MIFI OI
Suffix	Small Tip	Micro Fidelity & Extra Power	Open-Irrigated	Micro Fidelity & Open-Irrigated
Tip Length/Size	4mm/7F	8mm/8F	4mm/7F	4.5mm/7.5F
Shaft Size	7F	7F	7.5F	7.5F
Compatible Sheath	8F	8.5F	8F	8F
Curves	Standard, Large	Standard, Large, Asymmetric	Standard, Large, Asymmetric	Standard, Large, Asymmetric
RHYTHMIA Navigation Enabled	◆	◆	◆	◆
Open-Irrigated	-	-	◆	◆
Mini-Electrodes	-	◆	-	◆
Connection Cable	M004 RARC010	M004 RARC010	M004 RARC010	M004 RARC010



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi


Select from these topics

- Nav-Enabled Portfolio
  - Nav Technology
  - Mapping Creation
  - Mechanical Design
  - Blazer Platform
  - Construction - Conventional
  - Construction - OI
  - IntellaNav Spec Summary
  - IntellaNav Summary
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## INTELLANAV™ Specification Summary

**A Versatile Portfolio for Customized Ablation Strategy with RHYTHMIA HDx™**

- Tracking accuracy
- Map creation
- Ablation tip technology
  - Open-irrigation
  - Mini-electrodes
  - Tip lengths
- Familiarity of BLAZER™



INTELLAMAP ORION™  
MAPPING CATHETER

INTELLANAV MIFI™ OI  
ABLATION CATHETER

INTELLANAV™ OI  
ABLATION CATHETER

INTELLANAV MIFI™ XP  
ABLATION CATHETER

INTELLANAV™ ST  
ABLATION CATHETER

BETTER THAN  
**1mm** tracking  
ACCURACY  
IN RHYTHMIA™ MAPPING SYSTEM CASES.¹

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
  - Catheter Overview
  - Technical Specs
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Catheter Overview



***Unparalleled Clarity. Cool Performance. Confident Navigation.***

MiFi – Mini-electrode technology

Total Tip Cooling™ Design

Dual cooling chambers uniformly cool the entire tip internally

Optimized flow pattern actively washes entire tip electrode externally

Consistent cooling throughout RF delivery

BLAZER™ Platform Design

Dependable performance



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters


Open Irrigation

MiFi

Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
  - Catheter Overview
  - Technical Specs
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST

## Technical Specifications



### Platform

Built on the bidirectional BLAZER™ Catheter Platform

### Compatibility

Compatible with the RHYTHMIA™ Mapping System, MAESTRO 4000™ Cardiac Ablation System and METRIQ™ Pump

### Tip Design

Open-irrigated with 6 irrigation holes, dual internal cooling chambers and 3 mini-electrodes (4.5mm/7F)

### Shaft Size

7.5F

### Compatible Sheath

8F

### Curves

Standard, Large & Asymmetric

### Electrode Spacing

Quadripolar 2.5/2.5/2.5mm

### Shaft Length

110cm

### Proximal Shaft Material

Metallic braided Pebax™

### Distal Shaft Material

Fiber reinforced braided Pebax

### Electrodes

Pt/Ir 10%

### Maximum Wattage

50 W

### Generator Control Mode

Power



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
  - Catheter Overview
  - Technical Specs
- IntellaNav MiFi XP
- IntellaNav ST

## Catheter Overview



***Cool Performance. Confident Navigation.***

Total Tip Cooling™ Design

Dual cooling chambers uniformly cool the entire tip internally

Optimized flow pattern actively washes entire tip electrode externally

Consistent cooling throughout RF delivery

BLAZER™ Platform Design

Dependable performance



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi



Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
  - Catheter Overview
  - Technical Specs
- IntellaNav MiFi XP
- IntellaNav ST

## Technical Specifications



### Platform

Built on the bidirectional BLAZER™ Catheter Platform

### Compatibility

Compatible with the RHYTHMIA™ Mapping System, MAESTRO 4000™ Cardiac Ablation System and METRIQ™ Pump

### Tip Design

Open-irrigated with 6 irrigation holes and dual internal cooling chambers (4mm/7F)

### Shaft Size

7.5F

### Compatible Sheath

8F

### Curves

Standard, Large & Asymmetric

### Electrode Spacing

Quadripolar 2.5/2.5/2.5mm

### Shaft Length

110cm

### Proximal Shaft Material

Metallic braided Pebax™

### Distal Shaft Material

Pellethane™

### Electrodes

Pt/Ir 10%

### Maximum Wattage

50 W

### Generator Control Mode

Power



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
  - Catheter Overview
  - Technical Specs
- IntellaNav ST

## Catheter Overview



***Unparalleled Clarity. Performance. Accuracy.***

MiFi – Mini-Electrode Technology

Trusted BLAZER™ Platform

Torqueability (1:1 movements)

Synergy with RHYTHMIA™ High Definition Mapping

Bi-Directional Steering

Fine Micro-movements





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi



Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
  - Catheter Overview
  - Technical Specs
- IntellaNav ST

## Technical Specifications



### Platform

Built on the bidirectional BLAZER™ Catheter Platform

### Compatibility

Compatible with the RHYTHMIA™ Mapping System, MAESTRO 3000 & 4000™ Cardiac Ablation System

### Tip Design

Solid tip with 3 mini-electrodes (8mm/8F)

### Shaft Size

7F

### Compatible Sheath

8.5F

### Curves

Standard, Large & Asymmetric

### Electrode Spacing

Quadripolar 2.5/2.5/2.5mm

### Shaft Length

110cm

### Proximal Shaft Material

Metallic braided Pebax™

### Distal Shaft Material

Fiber reinforced braided Pebax

### Electrodes

Pt/Ir 10%

### Maximum Wattage

100 W

### Generator Control Mode

Temperature or Power



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST
  - Catheter Overview
  - Technical Specs

## Catheter Overview



***Familiarity. Performance. Accuracy.***

Historical performance meets HD mapping

Trusted BLAZER™ Platform

Torqueability (1:1 movements)

Pushability

Fine Micro-movements

In-plane Steering

Bi-Directional Curves

Precision of small ablation tip

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi



Select from these topics

- Nav-Enabled Portfolio
- IntellaNav MiFi OI
- IntellaNav OI
- IntellaNav MiFi XP
- IntellaNav ST
  - Catheter Overview
  - Technical Specs

## Technical Specifications



### Platform

Built on the bidirectional BLAZER™ Catheter Platform

### Compatibility

Compatible with the RHYTHMIA™ Mapping System, MAESTRO 3000 & 4000™ Cardiac Ablation System

### Tip Design

Solid tip (4mm/7F)

### Shaft Size

7F

### Compatible Sheath

8F

### Curves

Standard & Large

### Electrode Spacing

Quadripolar 2.5/2.5/2.5mm

### Shaft Length

110cm

### Proximal Shaft Material

Metallic braided Pebax™

### Distal Shaft Material

Fiber reinforced braided Pebax

### Electrodes

Pt/Ir 10%

### Maximum Wattage

50 W

### Generator Control Mode

Temperature or Power





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters


Open Irrigation

MiFi

Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System

**INTELLANAV™ OI**



**INTELLANAV MIFI™ OI**





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Cool Performance

### Total Tip Cooling™ Technology


Designed to reduce potential of char, coagulum and thrombus

### Dynamic external cooling

Proximally directed exit flow actively cools the entire tip electrode

### Active internal cooling

Dynamic saline flow fills two chambers to cool areas prone to hot spots from within



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters


Open Irrigation

MiFi

Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Total Tip Cooling™ Design



**Dual  
Cooling  
Chambers**

**Targeted  
Proximal  
Cooling**

**Exit Flow  
Directed  
Proximally**

**Central  
Cooling  
Lumen**

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi




Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System


## Same Internal Open-Irrigated Design

### INTELLANAV™ OI



Same design as BLAZER™ OI

### INTELLANAV MIFI™ OI



Same design as INTELLATIP MIFI™ OI





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi


Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Demonstrating Performance Equivalency: Lesions

### BLAZER™ OI VS. INTELLATIP MIFI™ OI RESULTS

(16 lesion samples per catheter per orientation)



### BLAZER™ OI & INTELLATIP MIFI™ OI

**Lesion dimensions (width & depth) are not clinically different.**

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics


- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Why open-irrigated design matters

Conventional open-irrigated catheters are designed to **reduce tip temperature** and enable delivery of higher power

But, tip cooling is **not uniform**, which can result in formation of **proximal hot spots**<sup>1</sup>

Cross Section of Conventional Open-Irrigated Catheter




Infrared Imaging

100.0°C



27.0°C



- Proximal heating
- Electrode
- Irrigation flow
- Lesion forming in tissue

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics


- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Total Tip Cooling™ Visualization


**Cross Section of Conventional Open-Irrigated Ablation Catheter**




**Infrared Imaging (20W, ~30 seconds)**




**External Flow Visualization (17 ml/min)**



**Cross Section of BLAZER™ Open-Irrigated Ablation Catheter**



**Infrared Imaging (20W, ~30 seconds)**



**External Flow Visualization (17 ml/min)**



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

● Cool Performance

- Total Tip Cooling™
- Internal OI Design
- Demonstrating Equivalency
- Why this matters
- Visualization
- Cooling Profile
- Technology
- Competitive Landscape
- Thermal Testing
- Temperature Sensing

○ Clinical Science

○ Open-Irrigated System


## Cooling Profile of BLAZER™ OPEN-IRRIGATED Catheter

### Design Goals

- Uniform cooling of tip electrode
- Cooling fluid actively washes tip electrode

### IR Thermography

(25W, 17 ml/min irrigation)



100.0°C

27.0°C



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
- Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System


## Total Tip Cooling™ Technology for Effective Power Delivery

### Efficient Tissue/Tip Interface Cooling

Promotes lesion creation while avoiding rapid tissue temperature rise

### BLAZER™ OI

Effective power delivery with steady tissue temperature rise



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Competitive Landscape

**ThermoCool™ 1st Generation vs. ThermoCool™ SF @30 seconds<sup>1</sup>**

**BLAZER™ OPEN-IRRIGATED vs. ThermoCool™ SF @30 seconds<sup>1</sup>**

Data based on bench testing conducted by Boston Scientific.


### Maximum Temperature during Ablation (20W Power Control)

#### ThermoCool SF

Bench testing showed rapid internal tissue temperature rise

#### ThermoCool

Bench testing showed power delivery with steady tissue temperature rise



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

● Cool Performance

○ Total Tip Cooling™

○ Internal OI Design

○ Demonstrating Equivalency

○ Why this matters

○ Visualization

○ Cooling Profile

○ Technology

● Competitive Landscape

○ Thermal Testing

○ Temperature Sensing

○ Clinical Science

○ Open-Irrigated System

## Competitive Landscape

ThermoCool™ 1st Generation vs. ThermoCool™ SF @30 seconds<sup>1</sup>

BLAZER™ OPEN-IRRIGATED vs. ThermoCool™ SF @30 seconds<sup>1</sup>

Data based on bench testing conducted by Boston Scientific.

### Maximum Temperature during Ablation (20W Power Control)

#### ThermoCool SF

Bench testing showed rapid internal tissue temperature rise

#### BLAZER OI

Bench testing showed power delivery with steady tissue temperature rise



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

● Cool Performance


- Total Tip Cooling™
- Internal OI Design
- Demonstrating Equivalency
- Why this matters
- Visualization
- Cooling Profile
- Technology
- Competitive Landscape
- Thermal Testing
- Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Thermal Testing | No compromise to cooling or power delivery

Data based on bench testing conducted by Boston Scientific.


### BLAZER™ OPEN-IRRIGATED Catheter vs. Competitive Catheters @30 seconds<sup>1</sup>

**BLAZER OPEN-IRRIGATED Catheter**




Saw total tip cooling with no hotspots  
Power delivery without rapid tissue temperature rise

**ThermoCool™ Catheter**



Tip cooling was not uniform, which can result in the formation of proximal hotspots  
Power delivery without rapid tissue temperature rise

**ThermoCool™ SF Low Flow Catheter**



Efficient tip cooling for no hot spots  
Saw rapid tissue temperature rise



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

## Select from these topics

- Cool Performance
  - Total Tip Cooling™
  - Internal OI Design
  - Demonstrating Equivalency
  - Why this matters
  - Visualization
  - Cooling Profile
  - Technology
  - Competitive Landscape
  - Thermal Testing
  - Temperature Sensing
- Clinical Science
- Open-Irrigated System

## Temperature Sensing

- INTELLANAV™ OI uses the same thermocouple and distal tip design as BLAZER™ OI
- The thermocouple is located in the distal tip and it is slightly offset with respect to the central axis of the catheter
- Physicians may notice a lower temperature during ablation
- Temperature reading is not indicative of effective ablation



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
  - Moreno et al.
  - Guerra et al.
  - Tip Temperature
  - Key Takeaways
- Open-Irrigated System

## Moreno et al. Thermodynamic Evaluation of 4 Open-Irrigated Catheters<sup>1</sup>

This study was an experimental assessment of four open-irrigated catheters comparing lesion size, safety and heat transfer. The thigh lesion model was performed on six anesthetized pigs.

### Thermodynamic comparison of four Open-Irrigated Catheters

- ThermoCool™
- ThermoCool™ SF
- BLAZER™ OPEN-IRRIGATED
- CoolFlex™

### Porcine thigh model

- 30W (60 sec) - assessing lesion morphology both perpendicularly and tangentially
- 50W (180 sec) – evaluating the propensity for deep-tissue overheating to compare “pop” rates
- Thermal assessment of lesion generation (20W, 60 sec)

### Results

- Newer catheters showed lower temp readings compared to ThermoCool
- ThermoCool SF pushed maximum thermal effect deeper, created largest lesions at 8mL/min (perpendicular)
- CoolFlex readily induced steam pops at 50W (perpendicular)
- No major differences found between BLAZER OI and ThermoCool (perpendicular)

LEARN MORE



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

Cool Performance

Clinical Science

Moreno et al.

Guerra et al.


Tip Temperature

Key Takeaways

Open-Irrigated System

## Moreno et al. Thermodynamic Evaluation of 4 Open-Irrigated Catheters<sup>1</sup>

**30W, 60 second lesions**



- No major differences were seen in maximum depth of lesions.
- CoolFlex showed a trend to lower volume lesions compared with BLAZER OI (P=0.09).
- No coagulum formation or thrombus occurred in any application.
- No differences seen in tangential (parallel) orientation (not shown).

← LEARN MORE →

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics

- Cool Performance
- Clinical Science
  - Moreno et al.
  - Guerra et al.
- Tip Temperature
- Key Takeaways
- Open-Irrigated System

## Moreno et al. Thermodynamic Evaluation of 4 Open-Irrigated Catheters<sup>1</sup>

Time-To-Pop Assessment (50 W, up to 180 seconds)						
	50W, 30 mL/min			50W, 15 mL/min	Global	
Perpendicular	BWThermoCool™	SJM CoolFlex™	BSC BLAZER™ OI	BW SF-30	BW SF-15	P Value
Tip Temperature (°C) Mean	38.4 ± 3.2	41.2 ± 6.5	34.8 ± 2.5	33 ± 1.4	33.4 ± 0.8	<0.001
Tip Temperature (°C) Maximum	43.9 ± 4.8	48.2 ± 14.6	42 ± 11.8	37.9 ± 1.8	36.8 ± 1.2	0.008
Pop before 30s	1 (8%)	8 (67%)	0 (0%)	2 (17%)	3 (25%)	
Pop before 60s	3 (25%)	12 (100%)	4 (33%)	2 (17%)	6 (50%)	
Pop before 90s	9 (75%)	12 (100%)	9 (75%)	8 (67%)	8 (67%)	
Pop before 180s	12 (100%)	12 (100%)	12 (100%)	9 (75%)	11 (91%)	





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
  - Moreno et al.
  - Guerra et al.
  - Tip Temperature
  - Key Takeaways
- Open-Irrigated System

## Guerra et al. In Vitro Comparison of 6 Open-Irrigated Catheters<sup>1</sup>

**The purpose of this study was to compare the lesion size and potential complications produced by commercially available open-irrigated catheters in an in vitro porcine heart model.**

- Catheters tested include: ThermoCool™, BLAZER™ OI, Therapy™ Cool Path™, Therapy™ Cool Path Duo™, ThermoCool™ SF and Therapy™ Cool Flex™

### **Total of 601 lesions were made in 26 in vitro preparations**

- 20 & 35 W, 30 & 60 sec
- Static flow rate of 13mL/min for all catheters

**Tip Temperature profile showed significant differences between the catheters (p < 0.001) with ThermoCool registering the lowest**

### **Lesion Volume**

- No significant differences between catheters found, regardless of power and flow conditions

### **Complications**

- Char mostly occurred at low-flow conditions
- Steam pops only occurred in the 35W, 60-second setting, being highest with Therapy Cool Path and ThermoCool SF

LEARN MORE



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

Cool Performance

Clinical Science

Moreno et al.

Guerra et al.


Tip Temperature

Key Takeaways

Open-Irrigated System

## Guerra et al. In Vitro Comparison of 6 Open-Irrigated Catheters<sup>1</sup>

Mean catheter tip temperature during the RF applications at the 2 tested powers and blood flow conditions. Bars represent: ± SEM.



At 20 W, when the blood flow increased a paradoxical response with an increase in tip temperature was seen with both ThermoCool™ SF and CoolFlex™.

- Because cooling effect is more homogenous, local fluid output per hole is lower, and, as a result, the irrigating saline flow is not able to counteract the higher blood flow.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters


Open Irrigation

MiFi

Select from these topics

- Cool Performance
- Clinical Science
  - Moreno et al.
  - Guerra et al.
  - Tip Temperature
  - Key Takeaways
- Open-Irrigated System

## Importance of Tip Temperature



Bench testing found, BLAZER™ OI had a consistently cool tip throughout RF delivery with steady tissue temperature rise



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



## Select from these topics

- Cool Performance
- Clinical Science
  - Moreno et al.
  - Guerra et al.
  - Tip Temperature
  - Key Takeaways
- Open-Irrigated System

## Key Takeaways:

### Moreno et al.<sup>1</sup>

- Catheter design has impacts on lesion creation and safety profile
- CoolFlex™ induced steam pops at higher powers
- ThermoCool™ SF pushed thermal effect deeper, creating larger lesions
- BLAZER™ OI did not induce any steam pops at 30 seconds

### Guerra et al.<sup>2</sup>

- Tip temperature profiles were statistically different
- No lesion creation difference seen among catheters
- Char occurred mostly in low flow conditions
- Newer “low flow” catheters saw an increase in temp in high blood flow
- BLAZER OI demonstrated a consistent cooling profile





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Introducing the MAESTRO 4000™ & METRIQ™ Pump

### Empowered Irrigated Ablation Procedures



Efficiency. Performance. Integration.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Open-Irrigated System

### Empower Efficiency


- Intelligent User Interface
- IntellaSight™ Infusion Monitoring
- Customizable Memory Buttons

### Empower Performance

- Bubble and Occlusion Detection
- Automatic Titration

### Empower Integration

- Cutting-Edge Catheter Portfolio
- Boston Scientific EP Lab



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Open-Irrigated System Components | MAESTRO 4000™ Cardiac Ablation System

### Seamless, automatic communication

- Flexible setup and easy operation
- Accurate, clear diagnostic message display
- Compatible with Boston Scientific's full portfolio of therapeutic and diagnostic catheters



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Open-Irrigated System Components | METRIQ™ Irrigation Pump

### Vital, real-time procedure information

- Large, easy-to-read display screen with intelligent user interface
- Intuitive, easy-to-interpret diagnostic messages
- Bubble and occlusion detection





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Intelligent User Interface

- Comprehensive, real-time diagnostic information on one screen
- Large, easy-to-read display can be viewed from a distance
- Quick, intuitive menu navigation



IntelliSight™ Infusion Monitoring

Customizable Memory Buttons



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## IntellaSight™ Infusion Monitoring

Provides real-time feedback on five different saline assessments

	Item	Description
1	Volume Remaining	Total fluid left in bag
2	Volume Infused	Total fluid dispensed while catheter is in patient
3	Volume Dispensed	Total fluid flow since start of procedure
4	New Saline Bag	Standard saline bag sizes from 500 mL to 2000 mL
5	Low Fluid Warning	Set alert message to display at 100 mL to 500 mL remaining

Automatic Mode		
Volume Remaining	1900 mL	2
Volume Infused	80 mL	
Volume Dispensed	100 mL	17
New Saline Bag?	YES	
<b>New Procedure?</b>	<b>YES</b>	30
Menu Page 1 of 3		

Automatic Mode		
<b>Saline Bag Size?</b>	<b>2000 mL</b>	2
Low Fluid Warning:	100 mL	
Fluid Vol. Display:	Remaining	17
Min. Temp. Drop:	2 deg. C	
To Trigger Hi Flow:	30 Watts	30
Menu Page 2 of 3		



IntellaSight™ Infusion Monitoring

Customizable Memory Buttons



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi




## Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Customizable Memory Buttons

- Programmable flow and menu settings can be quickly accessed for future procedures
- Store and retrieve procedure parameters for up to three different scenarios



### To store flow and menu settings

- Press and hold button for more than two seconds
- Pump will sound a short beep when button is pressed and a long beep when settings are stored
- LED will flash

### To retrieve flow and menu settings

- Press and then release button
- Pump sounds a short beep and the stored parameter will display

**White space can be used to label button**



IntellaSight™ Infusion Monitoring

Customizable Memory Buttons



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Bubble and Occlusion Detection

- Reliable sensor technology designed to prevent air infusion and occlusion
- Simple placement process to ensure accurate tubing alignment
- Smaller tubing designed to effectively clear bubbles and increase flow





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- Cool Performance
- Clinical Science
- Open-Irrigated System
  - Systems
  - Overview
  - MAESTRO 4000
  - METRIQ
  - User Interface
  - Detection
  - Titration

## Automatic Titration

- Intelligent, automatic titration ensures optimal power-to-fluid control
- Instant, clear display of titration status



Example Scenario: METRIQ™ Pump delivers saline at hi-flow (30 mL/min) when the MAESTRO 4000™ Controller power level is set to 30 Watts





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters


Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback

**INTELLANAV MIFI™ XP**



**INTELLANAV MIFI™ OI**





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
  - Unique Design
  - Unparalleled Clarity
  - Seeing Is Believing
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback

## Unique Catheter Design

### 3 Sophisticated Mini-Electrodes

- Enable localized recording of a small area
- Deliver signals of unparalleled clarity
- Allow multiple channels for highly localized EGMs



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi




Select from these topics

- MiFi Mini-Electrodes
  - Unique Design
  - Unparalleled Clarity
  - Seeing Is Believing
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback

## Unparalleled Clarity | The Truth with MicroFidelity

**Conventional Electrodes  
(bipolar mapping/pacing)**



**Mini-Electrodes  
(bipolar mapping/pacing)**



■ Center of Mapping / Pacing (bipolar)

■ Center of Ablation





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- MiFi Mini-Electrodes
  - Unique Design
  - Unparalleled Clarity
  - Seeing Is Believing
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback

## Map creation with INTELLANAV MIFI™ XP

### True Tip Location

MiFi sensitivity can show contact via EGMs

- Example 1 – concept animation
  - Example 2 – EGM strip
- MiFi can show when the catheter tip is:
- In the ventricle
  - In the atrium

Tip location in relation to a line of block

### True Tissue Assessment

- Viable or Non Viable Tissue
- vMap™ to assess block on the CTI
- Map double potentials
- Look for gaps
- Pacing thresholds to assess tissue

### True Ablation Feedback

- EGM attenuation
- Amplitude reduction greater than common bi-pole



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters



Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
  - Visual
  - Electrograms
  - All V on MIFI
  - A & V on MIFI
  - High Resolution Mapping
- True Tissue Assessment
- True Ablation Feedback

## True Tip Location



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi


Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
  - Visual
  - Electrograms
  - All V on MIFI
  - A & V on MIFI
  - High Resolution Mapping
- True Tissue Assessment
- True Ablation Feedback

## True Tip Location | Contact assessment with MiFi electrograms

### Distinguishing varying degrees of contact

- Tracing below shows respiratory movement leading to intermittent contact
- MiFis show not only “no contact” but also the transition from “good contact” to “no contact”
- Conventional bipoles (MAP-PROX and MAP-DIST) are dominated by far-field activity and suggest adequate contact throughout the recording.



**MiFi electrograms reveal more information than just contact or no contact.**

**They also reveal stability of contact.**



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
  - Visual
  - Electrograms
  - All V on MiFi
  - A & V on MiFi
  - High Resolution Mapping
- True Tissue Assessment
- True Ablation Feedback

True Tip Location | All V on MiFi

Mini-Electrode EGMs clearly demonstrate when the catheter tip is in the ventricle




Image courtesy of Kevin Makati, MD, St. Joseph's Hospital, Tampa, FL.





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
  - Visual
  - Electrograms
  - All V on MIFI
  - A & V on MIFI
  - High Resolution Mapping
- True Tissue Assessment
- True Ablation Feedback

True Tip Location | A & V on MiFi

Mini-Electrode EGMs clearly demonstrate when the catheter tip is entering the atrium




Image courtesy of Kevin Makati, MD, St. Joseph's Hospital, Tampa, FL.





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation



MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
  - Visual
  - Electrograms
  - All V on MIFI
  - A & V on MIFI
  - High Resolution Mapping
- True Tissue Assessment
- True Ablation Feedback

## Accurate High-Resolution Mapping

Maps generated with the mini-electrodes provide accurate information of tip location in relation to the line of block



Conventional electrode recordings suggest the tip electrode is located on the proximal side of the line of block (first potential of double potentials), but the mini-electrode recording shows that the tip electrode is actually located on the distal side (second potential of double potentials).

Nakagawa H, Harlev D, Koblish J, et al. Mini recording electrodes within a conventional 4.5 mm tip electrode improves ablation catheter mapping resolution. Poster session PO06-72 presented at 2012 Heart Rhythm Society, Boston, MA. (Right Atrial canine model, n=4. Operator utilized RHYTHMIA Mapping System with 4.5 mm Open-Irrigated IntellaTip MiFi catheter). Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
  - Tissue Substrate Identification
  - INTELLANAV MIFI & vMap
  - Mapping Double Potentials
  - Gap Identification
  - True Tissue Assessment
- True Ablation Feedback

## Tissue Substrate Identification

**INTELLANAV MIFI™ XP technology provided higher specificity and sensitivity in predicting atrial fibrosis and identifying abnormal substrate**




Image courtesy of Kevin Makati, MD, St. Joseph's Hospital, Tampa, FL.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi



Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
  - Tissue Substrate Identification
  - INTELLANAV MIFI & vMap
  - Mapping Double Potentials
  - Gap Identification
  - True Tissue Assessment
- True Ablation Feedback

## INTELLANAV MIFI™ XP & vMap™

**INTELLANAV MIFI XP provided clarity and accuracy to quickly vMap and assess bi-directional block**




Image courtesy of Matt Ostrom, MD, Torrance Memorial Hospital, Torrance, CA.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
  - Tissue Substrate Identification
  - INTELLANAV MIFI & vMap
  - Mapping Double Potentials
  - Gap Identification
  - True Tissue Assessment
- True Ablation Feedback

## Mapping Double Potentials Along the CTI Ablation Line




Image courtesy of Kevin Makati, MD, St. Joseph's Hospital, Tampa, FL.





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
  - Tissue Substrate Identification
  - INTELLANAV MIFI & vMap
  - Mapping Double Potentials
  - Gap Identification
  - True Tissue Assessment
- True Ablation Feedback

## Gap Identification

**High-resolution recording provided the sensitivity needed for quick and accurate detection of viable tissue**



1. Maddox W. The IntellaTip MiFi XP Ablation Catheter: Thoughts from a young electrophysiologist. Presented at 2013 Boston Scientific National Sales Meeting, Orlando FL. EP-222201-AA. Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.





# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi


Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
  - Tissue Substrate Identification
  - INTELLANAV MIFI & vMap
  - Mapping Double Potentials
  - Gap Identification
  - True Tissue Assessment
- True Ablation Feedback

## True Tissue Assessment

- Mini-electrodes exhibited marked threshold and EGM changes in the presence of transmural lesions in the atria
- The smaller surface area of the mini-electrodes compared to the 4.5-mm or 8-mm tip electrode yielded higher current density depolarizing a smaller tissue mass, leading to a lower pacing threshold

### ME Pacing can help assess viable vs non-viable tissue



Avital B, Horbal P, Vance D, et al. Maximal electrogram attenuation recorded from mini electrodes embedded on 4.5-mm irrigated and 8-mm nonirrigated catheters signifies lesion maturation. *J Cardiovasc Electrophysiol.* 2015 Feb;26(2):192-202. Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi




Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback
  - Clarity in EGM Attenuation
  - Lesion Maturation

## Clarity in EGM Attenuation

**Mini-electrodes demonstrated significant amplitude reduction and signal clarity during ablation as compared to unipolar and bipolar signals**



Case images courtesy of W. Jackman, MD, University of Oklahoma Health Sciences Center. Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.



# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation


MiFi

Select from these topics

- MiFi Mini-Electrodes
- True Tip Location
- True Tissue Assessment
- True Ablation Feedback
  - Clarity in EGM Attenuation
  - Lesion Maturation

## Lesion Maturation

- **EGM amplitude reduction, post ablation, was greater when measured with mini-electrodes compared to conventional 8 mm tip<sup>1</sup>**
- **EGM amplitude reduction on the mini-electrodes, post ablation, was correlated to transmuralit<sup>1</sup>**



1. Price A, Leshen Z, Hansen J, et al. Novel ablation catheter technology that improves mapping resolution and monitoring of lesion maturation. *J Innovat Cardiac Rhythm Manag.* 2012;3:599-609. (Canine model, n=7). Results from case studies are not necessarily predictive of results in other cases. Results in other cases may vary.

# NAV-ENABLED ABLATION CATHETERS

Overview

INTELLANAV™ Catheters

Open Irrigation

MiFi

7F = 2.33mm | 7.5F = 2.5 mm | 8F = 2.6mm | 8.5F = 2.33 mm

All cited trademarks are the property of their respective owners. CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. This material not intended for use in France.

