

NEUROMODULATION Radiofrequency Ablation

Product Catalogue

2019

Table of Contents

Generators	3
Generator Accessories	4
Nitinol Reusable Electrodes	5
Stainless Steel Reusable Electrodes	6
Disposable Electrodes	7
Unified Disposable Electrode/Cannula	8
Disposable Injection Electrodes and Needles	10
Disposable Procedure Kits	11
Straight Sharp Cannulae	12
Curved Sharp Cannulae	14
Curved Blunt Cannulae	15
Sidekick™ Cannulae	16
Halyard Curved Cannula	16
Disposable Ground Pads	17
Disposable Palisade™ Guideblocks and Kits	18
Specialty Electrodes and Cannulae Kits	19
Connection Cables and Adapters	20




Generators

PRODUCT NUMBER	DESCRIPTION
RFG-4-220V	<p>G4™ RF Generator: 220 V Kit</p> <p><i>Package Contents</i></p> <ul style="list-style-type: none"> (1) G4 Generator (1) Medical Grade Power Cord (1) RFG-TP Generator Output Test Plug (1) CB112-TC Connection Cable for Reusable Electrodes
RFG-4-240V	<p>G4™ RF Generator: 240 V Kit</p> <p><i>Package Contents</i></p> <ul style="list-style-type: none"> (1) G4 Generator (1) Medical Grade Power Cord (1) RFG-TP Generator Output Test Plug (1) CB112-TC Connection Cable for Reusable Electrodes





Generator Accessories

PRODUCT NUMBER	DESCRIPTION
COR-10A-NA	<p>Medical Grade Power Cord for 120V generator For use in: USA, Canada, Columbia, Costa Rica, Guatemala, Honduras, Mexico, Brazil, Colombia, Panama, Taiwan</p> 
COR-10A-IN	<p>Medical Grade Power Cord for 220V generator For use in: Denmark, France, Israel, Korea, Argentina, Austria, Belgium, Brazil, Bulgaria, Germany, Hong Kong, Luxembourg, Netherlands, Russia, Spain, Turkey, Ukraine, Vietnam</p> <p>Medical Grade Power Cord for 240V generator For use in: Iraq, Pakistan, Poland, Bahrain, India</p>
COR-10A-IT	<p>Medical Grade Power Cord for 220V generator For use in: Italy</p>
COR-10A-UK	<p>Medical Grade Power Cord for 220V generator For use in: UK, Dubai, Saudi Arabia, South Africa</p> <p>Medical Grade Power Cord for 240V generator For use in: Singapore, UK, Ireland</p>
COR-10A-SW	<p>Medical Grade Power Cord for 220V generator For use in: Switzerland</p>
FS-1	<p>Foot Switch</p> 
RFG-TP	<p>Generator Output Test Plug</p> 
SC-1	<p>Storage and Carrying Case for RFG-1A or RFG-1B</p>
SC-2	<p>Storage and Carrying Case for G4</p>



Nitinol Reusable Electrode

PRODUCT NUMBER	DESCRIPTION
	Nitinol Reusable Electrode Please note: These electrodes require connection cable CB112-TC 
TCN-10	Nitinol TC Electrode: 10 cm
TCN-10-3M	Nitinol TC Electrode: 10 cm, 3m Cable
TCN-15	Nitinol TC Electrode: 15 cm
TCN-15-3M	Nitinol TC Electrode: 15 cm, 3m Cable
TCN-20	Nitinol TC Electrode: 20 cm
TCN-20-3M	Nitinol TC Electrode: 20 cm, 3m Cable
TCN-5	Nitinol TC Electrode: 5 cm
TCN-5-3M	Nitinol TC Electrode: 5 cm, 3m Cable
TCN-CASE	Sterilizing Case for CSK Electrode

Stainless Steel Reusable Electrodes

PRODUCT NUMBER	DESCRIPTION
	Stainless Steel Reusable Electrode Please note: These electrodes require connection cable CB112-TC
CSK-TC5	Stainless Steel Reusable Electrode – 5 cm
CSK-TC10	Stainless Steel Reusable Electrode – 10 cm
CSK-TC15	Stainless Steel Reusable Electrode – 15 cm
CSK-TC20	Stainless Steel Reusable Electrode – 20 cm
CB112-TC	Connection Cable for Reusable Electrodes 
	Stainless Steel Reusable Electrode with integrated 3m Cable 
CSK-TC5-3M	Stainless Steel Reusable Electrode – 5 cm with integrated 3m cable
CSK-TC10-3M	Stainless Steel Reusable Electrode – 10 cm with integrated 3m cable
CSK-TC15-3M	Stainless Steel Reusable Electrode – 15 cm with integrated 3m cable
CSK-CASE	Sterilizing Case for CSK Electrode 

Disposable Electrodes




PRODUCT NUMBER	DESCRIPTION
	<p>Disposable Electrode</p> <p><i>Package Contents</i> (10) Disposable Electrodes</p> <p>Please note: These electrodes require connection cable CB114-TC</p> 
TCD-5-P	Disposable Electrode – 5 cm
TCD-10-P	Disposable Electrode – 10 cm
TCD-15-P	Disposable Electrode – 15 cm
TCD-20-P	Disposable Electrode – 20 cm
	<p>Disposable Electrode with integrated 3m Cable</p> <p><i>Package Contents</i> (10) Disposable Electrodes</p> <p>Please note: These electrodes require adapter CB114-ADP or connection cable CB114-TC</p> 
TCD-5-3M-P	Disposable Electrode – 5 cm with integrated 3m cable
TCD-10-3M-P	Disposable Electrode – 10 cm with integrated 3m cable
TCD-15-3M-P	Disposable Electrode – 15 cm with integrated 3m cable
CB114-ADP	Adapter Plug for Disposable Electrodes
CB114-TC	Connection Cable for Disposable Electrodes

Unified Disposable Electrode/Cannula




PRODUCT NUMBER	DESCRIPTION
	Unified Disposable Electrode/Cannula <i>Package Contents</i> (10) Unified Disposable Electrode/Cannula Please note: All CU/CUR electrodes require connection cable CB114-TC
	Unified Disposable Electrode/Cannula Straight Tip
CU-6522-P	Length: 6 cm – Tip: 5 mm – Gauge: 22 ga
CU-61022-P	Length: 6 cm – Tip: 10 mm – Gauge: 22 ga
CU-10522-P	Length: 10 cm – Tip: 5 mm – Gauge: 22 ga
CU-101018-P	Length: 10 cm – Tip: 10 mm – Gauge: 18 ga
CU-101022-P	Length: 10 cm – Tip: 10 mm – Gauge: 22 ga
CU-15522-P	Length: 15 cm – Tip: 5 mm – Gauge: 22 ga
	Unified Disposable Electrode/Cannula Straight Tip – Echo RF™
CU-6520U-P	Length: 6 cm – Tip: 5 mm – Gauge: 20 ga
CU-10520U-P	Length: 10 cm – Tip: 5 mm – Gauge: 20 ga
CU-101020U-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
CU-15520U-P	Length: 15 cm – Tip: 5 mm – Gauge: 20 ga
CU-151020U-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga






Unified Disposable Electrode/Cannula

PRODUCT NUMBER	DESCRIPTION
	Unified Disposable Electrode/Cannula Curved Tip 
CUR-6518-P	Length: 6 cm – Tip: 5 mm – Gauge: 18 ga
CUR-101018-P	Length: 10 cm – Tip: 10 mm – Gauge: 18 ga
CUR-151018-P	Length: 15 cm – Tip: 10 mm – Gauge: 18 ga
	Unified Disposable Electrode/Cannula Curved Tip – Echo RF™ 
CUR-101020U-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
CUR-151020U-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga
CB114-TC	Connection Cable for Disposable Electrodes 



Disposable Injection Electrodes and Needles

PRODUCT NUMBER	DESCRIPTION
	<p>Voltage Controlled Injection Electrode (Non-temperature Controlled)</p> <p><i>Package Contents</i> (10) Voltage Controlled Injection Electrodes</p> <p>Please note: These electrodes require connection cable CB116 or CB116-4</p> 
CR-6-P	Length: 6 cm – Tip: 5 mm
CR-6522U-P	Length: 6 cm – Tip: 5 mm – Echo RF
CR-10-P	Length: 10 cm – Tip: 5 mm
CR-10522U-P	Length: 10 cm – Tip: 5 mm – Echo RF
CR-15-P	Length: 15 cm – Tip: 5 mm
	<p>Stimulation Guided Nerve Block Needle</p> <p><i>Package Contents</i> (10) Stimulation Guided Nerve Block Needles</p> <p>Please note: These needles require connection cable CB116 or CB116-4</p> 
CP-6-P	Length: 6 cm – Tip: 1 mm
CP-10-P	Length: 10 cm – Tip: 1 mm
	<p>Disposable Nerve Block Needle</p> <p><i>Package Contents</i> (10) Disposable Nerve Block Needles</p> 
CN-6-P	Length: 6 cm
CN-10-P	Length: 10 cm
CN-15-P	Length: 15 cm

Disposable Injection Electrodes and Needles

PRODUCT NUMBER	DESCRIPTION	
	Disposable RF Injection Electrode Connection Cables	
CB116	Single Injection Electrode Cable	
CB116-4	Four Injection Electrode Cable	
CB103-1A	Injection Electrode to Alligator Clip Adapter	

Disposable Procedure Kits

PRODUCT NUMBER	DESCRIPTION	
BDK-15	<p>Bipolar Disc Kit</p> <p><i>Package Contents</i></p> <p>(2) CC152020 Straight Sharp Cannulae</p> <p>(2) TCD-15 15cm Disposable Electrodes</p> <p>Please note: These electrodes require connection cable CB114-TC</p>	
CB114-TC	Connection Cable for Disposable Electrodes	

Straight Sharp Cannulae

PRODUCT NUMBER	DESCRIPTION
	Straight Sharp Cannula <i>Package Contents</i> (10) Straight Sharp Cannulae
	
CC5222-P	Length: 5 cm – Tip: 2 mm – Gauge: 22 ga
CC5421-P	Length: 5 cm – Tip: 4 mm – Gauge: 21 ga
CC5422-P	Length: 5 cm – Tip: 4 mm – Gauge: 22 ga
CC51021-P	Length: 5 cm – Tip: 10 mm – Gauge: 21 ga
CC10222-P	Length: 10 cm – Tip: 2 mm – Gauge: 22 ga
CC10518-P	Length: 10 cm – Tip: 5 mm – Gauge: 18 ga
CC10520-P	Length: 10 cm – Tip: 5 mm – Gauge: 20 ga
CC10521-P	Length: 10 cm – Tip: 5 mm – Gauge: 21 ga
CC10522-P	Length: 10 cm – Tip: 5 mm – Gauge: 22 ga
CC101018-P	Length: 10 cm – Tip: 10 mm – Gauge: 18 ga
CC101020-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
CC101021-P	Length: 10 cm – Tip: 10 mm – Gauge: 21 ga
CC101022-P	Length: 10 cm – Tip: 10 mm – Gauge: 22 ga
CC15220-P	Length: 15 cm – Tip: 2 mm – Gauge: 20 ga
CC15518-P	Length: 15 cm – Tip: 5 mm – Gauge: 18 ga
CC15520-P	Length: 15 cm – Tip: 5 mm – Gauge: 20 ga
CC151018-P	Length: 15 cm – Tip: 10 mm – Gauge: 18 ga
CC151020-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga
CC151520-P	Length: 15 cm – Tip: 15 mm – Gauge: 20 ga
CC152020-P	Length: 15 cm – Tip: 20 mm – Gauge: 20 ga
CC20520-P	Length: 20 cm – Tip: 5 mm – Gauge: 20 ga
CC201020-P	Length: 20 cm – Tip: 10 mm – Gauge: 20 ga



Straight Sharp Cannulae


PRODUCT NUMBER	DESCRIPTION
	Straight Sharp Cannula – Echo RF™ <i>Package Contents</i> (10) Straight Sharp Cannulae – Echo RF
	
CC5521U-P / CC5521VU-P	Length: 5 cm – Tip: 5 mm – Gauge: 21 ga
CC51021U-P	Length: 5 cm – Tip: 10 mm – Gauge: 21 ga
CC10520U-P	Length: 10 cm – Tip: 5 mm – Gauge: 20 ga
CC10521U-P	Length: 10 cm – Tip: 5 mm – Gauge: 21 ga
CC101020U-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
CC101021U-P / CC101021VU-P	Length: 10 cm – Tip: 10 mm – Gauge: 21 ga
CC15520U-P	Length: 15 cm – Tip: 5 mm – Gauge: 20 ga
CC151020U-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga
CC151021VU-P	Length: 15 cm – Tip: 10 mm – Gauge: 21 ga

Straight Sharp Cannulae – Premium Siliconized

PRODUCT NUMBER	DESCRIPTION
	Straight Sharp Cannula – Premium Siliconized <i>Package Contents</i> (10) Straight Sharp Cannulae
CC5422Z-P	Length: 5 cm – Tip: 4 mm – Gauge: 22 ga – Premium Siliconized
CC10522Z-P	Length: 10 cm – Tip: 5 mm – Gauge: 22 ga – Premium Siliconized
CC101020Z-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga – Premium Siliconized
CC151020Z-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga – Premium Siliconized
CC101022Z-P	Length: 10 cm – Tip: 10 mm – Gauge: 22 ga – Premium Siliconized



Curved Sharp Cannulae

PRODUCT NUMBER	DESCRIPTION
	Curved Sharp Cannula <i>Package Contents</i> (10) Curved Sharp Cannulae
	
RFK-C5422S-P	Length: 5 cm – Tip: 4 mm – Gauge: 22 ga
RFK-C5616S-P	Length: 5 cm – Tip: 6 mm – Gauge: 16 ga
RFK-C51016S-P	Length: 5 cm – Tip: 10 mm – Gauge: 16 ga
RFK-C10518S-P	Length: 10 cm – Tip: 5 mm – Gauge: 18 ga
RFK-C10520S-P	Length: 10 cm – Tip: 5 mm – Gauge: 20 ga
RFK-C10522S-P	Length: 10 cm – Tip: 5 mm – Gauge: 22 ga
RFK-C10616S-P	Length: 10 cm – Tip: 6 mm – Gauge: 16 ga
RFK-C101016S-P	Length: 10 cm – Tip: 10 mm – Gauge: 16 ga
RFK-C101018S-P	Length: 10 cm – Tip: 10 mm – Gauge: 18 ga
RFK-C101020S-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
RFK-C15518S-P	Length: 15 cm – Tip: 5 mm – Gauge: 18 ga
RFK-C15520S-P	Length: 15 cm – Tip: 5 mm – Gauge: 20 ga
RFK-C15616S-P	Length: 15 cm – Tip: 6 mm – Gauge: 16 ga
RFK-C151016S-P	Length: 15 cm – Tip: 10 mm – Gauge: 16 ga
RFK-C151018S-P	Length: 15 cm – Tip: 10 mm – Gauge: 18 ga
RFK-C151020S-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga
RFK-C151520S-P	Length: 15 cm – Tip: 15 mm – Gauge: 20 ga
RFK-C20520S-P	Length: 20 cm – Tip: 5 mm – Gauge: 20 ga



Curved Sharp Cannulae – Premium Siliconized


PRODUCT NUMBER	DESCRIPTION
	Curved Sharp Cannula Premium Siliconized <i>Package Contents</i> (10) Curved Blunt Cannulae
RFK-C10522SZ-P	Length: 10 cm – Tip: 5 mm – Gauge: 22ga – Premium Siliconized
RFK-C101016VZ-P	Length: 10 cm – Tip: 10 mm – Gauge: 16ga – Premium Siliconized
RFK-C101018VZ-P	Length: 10 cm – Tip: 10 mm – Gauge: 18ga – Premium Siliconized
RFK-C101020VZ-P	Length: 10 cm – Tip: 10 mm – Gauge: 20ga – Premium Siliconized
RFK-C151018VZ-P	Length: 15 cm – Tip: 10 mm – Gauge: 18ga – Premium Siliconized
RFK-C151020VZ-P	Length: 15 cm – Tip: 10 mm – Gauge: 20ga – Premium Siliconized

Curved Blunt Cannulae


PRODUCT NUMBER	DESCRIPTION
	Curved Blunt Cannula <i>Package Contents</i> (10) Curved Blunt Cannulae
RFK-C10520B-P	Length: 10 cm – Tip: 5 mm – Gauge: 20 ga
RFK-C101020B-P	Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
RFK-C15520B-P	Length: 15 cm – Tip: 5 mm – Gauge: 20 ga
RFK-C151020B-P	Length: 15 cm – Tip: 10 mm – Gauge: 20 ga









Sidekick™ Cannulae

PRODUCT NUMBER	DESCRIPTION
	Sidekick Cannula <i>Package Contents</i> (10) Sidekick Cannulae
	
RFK-C101018S-ZK-P	Length: 10 cm – Tip: 10 mm – Gauge: 18ga
RFK-C101020S-ZK-P	Length: 10 cm – Tip: 10 mm – Gauge: 20ga
RFK-C151018S-ZK-P	Length: 15 cm – Tip: 10 mm – Gauge: 18ga
RFK-C151020S-ZK-P	Length: 15 cm – Tip: 10 mm – Gauge: 20ga




Halyard Curved Cannulae

PRODUCT NUMBER	DESCRIPTION
	Halyard Curved Cannula <i>Package Contents</i> (10) Halyard Curved Cannulae
	
	<i>Please note: These cannulae are designed to work with Halyard Curved electrodes only</i>
RFK21-54-5CS-P	Length: 54 mm – Tip: 5 mm – Gauge: 21 ga
RFK18-100-10CS-P	Length: 100 mm – Tip: 10 mm – Gauge: 18 ga
RFK20-100-10CS-P	Length: 100 mm – Tip: 10 mm – Gauge: 20 ga
RFK20-145-10CS-P	Length: 145 mm – Tip: 10 mm – Gauge: 20 ga





Disposable Ground Pads

PRODUCT NUMBER	DESCRIPTION	
4777	Disposable Ground pad w/cord <i>Package Contents</i> (1) Ground pad	
4777-5	Disposable Ground pad w/cord <i>Package Contents</i> (5) Ground pads	
2407	Disposable Ground pad cordless <i>Package Contents</i> (1) Ground pad	
2407-5	Disposable Ground pad cordless <i>Package Contents</i> (5) Ground pads	
CB105	Connection Cable for Cordless Ground Pads	
CB103-B	Ground Pad Cable for Non-Disposable Ground Plate	

Disposable Palisade™ Guideblocks and Kits

PRODUCT NUMBER	DESCRIPTION
	Disposable Palisade Guideblock <i>Package Contents</i> (10) Disposable Palisade Guideblocks 
PG-1018-P	Disposable Palisade Guideblock for 18 ga Cannulae
PG-1020-P	Disposable Palisade Guideblock for 20 ga Cannulae
	Disposable Palisade Guideblock Kit <i>Package Contents</i> (6) Straight Sharp Cannulae (4) TCD-10 Disposable Electrodes These electrodes require connection cable CB114-TC Disposable Palisade Guideblock to be ordered separately 
PALISADE-18	Included Cannulae: CC101018 - Length: 10 cm – Tip: 10 mm – Gauge: 18 ga
PALISADE-20	Included Cannulae: CC101020 - Length: 10 cm – Tip: 10 mm – Gauge: 20 ga
CB114-TC	Connection Cable for Disposable Electrodes 


















Specialty Electrodes and Cannulae Kits

PRODUCT NUMBER	DESCRIPTION
FLEX-K151320-5P	<p data-bbox="392 477 676 510">Flextrode™ Introducer Kit</p> <p data-bbox="392 544 751 712"><i>Package Contents</i> (5) Flextrode Disposable Electrodes Tip: 13 mm – Gauge: 20 ga (5) Flextrode Cannulae Tip: 15 cm – Gauge: 17 ga</p>  <p data-bbox="392 750 1145 779"><i>Please note: The Flextrode Introducer Kit requires connection cable CB114-TC</i></p>
	<p data-bbox="392 808 943 842">Reig-Cosman Epidural (RCE) Pulsed RF Electrode</p> <p data-bbox="392 875 711 1010"><i>Package Contents</i> (10) RCE Disposable Electrodes Length: 40 cm Length: 66 cm</p>  <p data-bbox="392 1043 1114 1072"><i>Please note: The Pulsed Electrode Kit requires connection cable CB114-TC</i></p>
RCE-E401519-P	<p data-bbox="392 1111 655 1144">(10) RCE 40 cm Electrodes</p> 
RCE-E661519-P	<p data-bbox="392 1184 655 1218">(10) RCE 66 cm Electrodes</p>
RCE-C916S-P	<p data-bbox="392 1258 691 1292">(10) RCE Introducer Cannulae</p>
CB114-TC	<p data-bbox="392 1332 890 1366">Connection Cable for Disposable Electrodes</p>  <p data-bbox="392 1400 1066 1429"><i>Please note: The Bipolar Disc Kit requires connection cable CB114-TC</i></p>

Connection Cables and Adapters

PRODUCT NUMBER	DESCRIPTION	
	Connection Cables and Adapters for Boston Scientific (Cosman) Generators	
CB112-TC	Connection Cable for Reusable Electrodes	
CB114-TC	Connection Cable for Disposable Electrodes	
CB114-ADP	Adapter Plug for Disposable Electrodes	
CB105	Connection Cable for Cordless Ground Pads	
CB116	Single Injection Electrode Cable	
CB116-4	Four Injection Electrode Cable	
CB103-1A	Injection Electrode to Alligator Clip Adapter	
CB103-B	Ground Pad Cable for Non-Disposable Ground Plate	
CB103-R	Alligator Clip Connection Cable	
CB113	RFG-1A/1B Bipolar Adapter Cable for Ground Pad Jack (Reusable Electrodes)	
CB130-TC	Neurotherm Simplicity Electrode to Cosman Generator Cable	
CB135-TC	Neurotherm Reusable Electrode to Cosman Generator Cable	
CB-SN	Smith & Nephew Spinecath Electrode to Cosman Generator Cable	

Connection Cables and Adapters

PRODUCT NUMBER	DESCRIPTION	
	Connection Cables and Adapters for Major Manufacturers Generators <i>Please note: Some adapters may require an additional connection cable for use with some generators</i>	
C125-TC	Adapter for Stryker Generators	
C120-TC	Adapter Cable for Halyard Single-Output Generators (Reusable Electrodes)	
C121-TC	Adapter Cable for Halyard Single-Output Generators (Disposable Electrodes)	
C124-TC	Adapter Cable for Halyard Four-Output Generators (Reusable Electrodes)	
C134-TC	Adapter Cable for Halyard Four-Output Generators (Disposable Electrodes)	
C126-TC	Adapter Cable for Smith & Nephew Generators (Reusable Electrodes)	
C127-TC	Adapter Cable for Smith & Nephew Generators (Disposable Electrodes)	
C115-TC	Adapter Cable for Neurotherm Generators (Reusable Electrodes)	
C123-TC	Adapter Cable for Neurotherm Generators (Disposable Electrodes)	
C110	Adapter Cable for Neurotherm Generators (Injection Electrodes)	
C112-TC	Adapter Cable for Radionics Generators (Reusable Electrodes)	
C114-TC	Adapter Cable for Radionics Generators (Disposable Electrodes)	
C116	Adapter Cable for Radionics Generators (Injection Electrodes)	
C116-4	Adapter Cable for Radionics Generators (Injection Electrodes x 4)	
C105	Ground Pad Adapter Cable for Halyard/Radionics Generators (DGP-PMC)	
C103-RTM	Radionics Thermistor Cable to Alligator Clip	
C104-TM	Radionics Thermistor Cable for Thermistor Electrode	



Certain Terms and Conditions

WARRANTY

- Radiofrequency Generators – Limited warranty for a period of 1 year from the date of purchase.
- Please see Limited Warranty Manual (91006760-01) for more detailed information.

RADIOFREQUENCY ABLATION CUSTOMER SERVICE

Please contact your local or direct customer service team by email or phone for more information.

TERMS

Terms are Net 30 Days. All shipments are F.O.B. from our distribution centers. Prices are subject to change without notice.

RETURNED GOODS

Permission to return merchandise is required for all returns and may be withheld at the sole discretion of Boston Scientific Neuromodulation. This authorization to return goods must be requested by calling or writing to your local/direct Boston Scientific Customer Service Team.

MERCHANDISE THAT WILL NOT BE ACCEPTED FOR RETURN:

- Expired product
- Merchandise in unsaleable condition (opened, damaged, tampered, etc.)
- Merchandise not listed in this Catalog, Price List, or any Supplement
- Merchandise sterilized or used by the customer
- Merchandise that is not packaged properly for shipping when returned
- Incomplete kits or boxes (unless defective)
- Unauthorized returns

SPECIFICATION CHANGES

As we are continually introducing products and improving our existing products, this Catalog does not purport to be complete, and the equipment supplied may not agree in all details with the description and specifications shown in the catalog. All product offerings, prices, and specifications are subject to change without notice.

NOTICE

Indications, warnings, contraindications, and instructions for use, as well as product warranty information, can be found in the product labeling supplied with each device. **EXCEPT AS SET FORTH IN THE PRODUCT LABELING SUPPLIED WITH THE DEVICE, BOSTON SCIENTIFIC NEUROMODULATION MAKES NO WARRANTIES, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, WITH RESPECT TO THE DEVICE.**

WARNING RE: SINGLE-USE DEVICES

Products marked as “single-use” are for single use only. Do not reuse, reprocess, or resterilize single-use products. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing, or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient. Boston Scientific Neuromodulation assumes no liability with respect to instruments reused, reprocessed, or resterilized, and makes no warranties, express or implied, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose, with respect to any such instrument.





Information contained herein is for distribution outside the U.S. only. Illustrations for information purposes – not indicative of actual size or clinical outcome.

All cited trademarks are the property of their respective owners. **CAUTION:** The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labeling supplied with each device. Information for the use only in countries with applicable health authority product registrations. Material not intended for use in France.

NM-539001-AB Printed in Germany by medicalvision.

**Boston
Scientific**
Advancing science for life™

www.bostonscientific.eu

© 2019 Boston Scientific Corporation
or its affiliates. All rights reserved.