

# ELI 150c

MORTARA'S NEXT GENERATION 12-LEAD ELECTROCARDIOGRAPHS

## PRODUCT FEATURES

- **Portable Solution**

Compact and lightweight, the ELI™ 150c electrocardiograph provides comprehensive functionality in an ultra-portable device.

- **High-Resolution Color LCD**

High-resolution color display provides real-time preview of 12-lead ECG and post-acquisition review of acquired ECG.

- **Distinguished Mortara VERITAS™ Resting ECG Interpretation Algorithm**

Widely recognized Mortara resting ECG interpretation algorithm uses gender specific and adult and pediatric criteria to provide a silent second opinion for resting ECG interpretation.

- **Choice of Wireless or Traditional ECG Acquisition**

ELI 150c offers a choice of either the innovative WAM™ wireless acquisition module or the AM12™ acquisition module. Both include replaceable lead wires, lead fail indicator, and remote control with buttons for ECG acquisition and rhythm printing.

- **Best 10**

ELI 150c automatically selects and displays the Best 10 seconds of ECG signal from the five minute segment of internal full disclosure, reducing clinical review time and eliminating the need for repeat ECGs.

- **Full Keyboard**

Alphanumeric elastomer keyboard features dedicated “one-touch” buttons for ECG acquisition, rhythm printing, and ECG transmission/order retrieval.


- **Information Exchange**

Experience the benefits of bidirectional communication via internal modem, LAN, wireless LAN, or GPRS mobile with Mortara's E-Scribe™ and Athena products, as well as third party EMR, HIS, and PACS systems via XML, PDF, UNIPRO32, and DICOM®32.



12-LEAD resting electrocardiograph

 **Mortara**





# ELI™ 150c ELECTROCARDIOGRAPH

## MORTARA INSTRUMENT, INC.

7865 North 86th Street  
Milwaukee, WI 53224  
U.S.A.

Tel: 414.354.1600  
Tel: 800.231.7437  
Service: 888.MORTARA  
Fax: 414.354.4760

## MORTARA RANGONI EUROPE, SRL

(European Headquarters)  
Via Cimarosa 103/105  
40033 Casalecchio di Reno (BO)  
Italy

Tel: +39.051.298.7811  
Fax: +39.051.613.3582

## MORTARA INSTRUMENT AUSTRALIA

PO Box 7568  
Baulkham Hills NSW 2153  
Unit 28, 9 Hoyle Avenue  
Castle Hill NSW 2154  
Australia

Tel: +61 2 8070 9303  
Fax: +61 2 9899 9478

[www.mortara.com](http://www.mortara.com)

ISO 13485 CERTIFIED



## WARRANTY + SERVICE

Mortara Instrument is committed to the highest level of customer support. Please contact us for the program which best suits your needs.

FEATURE	SPECIFICATION*
<b>INSTRUMENT TYPE</b>	12-lead electrocardiograph
<b>INPUT CHANNELS</b>	Simultaneous acquisition of all 12 leads
<b>STANDARD LEADS ACQUIRED</b>	I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6
<b>WAVEFORM DISPLAY</b>	Backlit, 1/4 VGA color LCD (320 x 240); 3, 4+4, or 6+6 lead presentation
<b>INPUT IMPEDANCE INPUT DYNAMIC RANGE ELECTRODE OFFSET TOLERANCE COMMON MODE REJECTION</b>	Meets or exceeds the requirements of ANSI/AAMI EC11
<b>PATIENT LEAKAGE CURRENT CHASSIS LEAKAGE CURRENT</b>	Meets or exceeds the requirements of ANSI/AAMI ES1
<b>DIGITAL SAMPLING RATE</b>	40,000 s/sec/channel used for pacemaker spike detection; 1000 s/sec/channel used for recording and analysis
<b>OPTIONAL FUNCTIONS</b>	Mortara VERITAS resting ECG interpretation algorithm with age and gender specific criteria; connectivity with bidirectional communication
<b>PAPER</b>	Perforated double Z-fold thermal paper; 108 mm (4") wide, 200 sheets
<b>THERMAL PRINTER</b>	Computer-controlled dot array; 8 dots/mm
<b>THERMAL PRINTER SPEEDS</b>	5, 10, 25, or 50 mm/s
<b>GAIN SETTINGS</b>	5, 10, or 20 mm/mV
<b>REPORT PRINT FORMATS</b>	Standard or Cabrera; 3, 3+1, 3+3, or 6 channel
<b>RHYTHM PRINT FORMATS</b>	3 or 6 channel with configurable lead groups
<b>KEYBOARD</b>	Elastomer keyboard with complete alphanumeric keys, soft-key menu, and dedicated function keys
<b>FREQUENCY RESPONSE</b>	0.05 to 300 Hz
<b>FILTERS</b>	High-performance baseline filter; AC interference filter 50/60 Hz; low-pass filters 40 Hz, 150 Hz, or 300 Hz
<b>A/D CONVERSION</b>	20 bits (1.17 microvolt LSB)
<b>DEVICE CLASSIFICATION</b>	Class I, Type CF defibrillation-proof applied parts
<b>ECG STORAGE</b>	Internal storage up to 100 ECGs; optional expanded up to 200 ECGs
<b>WEIGHT</b>	7.2 lbs. (3.3 kg) including battery (without paper)
<b>DIMENSIONS</b>	11.25 x 11.5 x 3.75" (29.2 x 30.5 x 10.2 cm)
<b>POWER REQUIREMENTS</b>	Universal AC power supply (100-240 VAC at 50/60 Hz) 110 VA; internal rechargeable battery

\*Specifications subject to change without notice.

AM12™, ELI™, E-Scribe™, VERITAS™ and WAM™ are trademarks of Mortara Instrument, Inc.

DICOM® is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.